




tepecanada.com
888-789-9928

PREVENTING PERI-IMPLANT COMPLICATIONS

Caring for implants



There are various types of implants which require proper **home care** maintenance. TePe offers a wide range of products to make cleaning easier for patients.





Special brushes

Patients may need a special brush to be able to clean the area where the implants meet the soft tissue.

TePe Implant Care™ has an **angled neck** which reaches the implant surfaces from the lingual. **TePe Implant Brush** has a slim brush head for easy cleaning of the buccal implant surfaces.



TePe Interdental Brush original.



Single tufted brush

A single tufted brush is ideal for precision cleaning of areas where it is difficult to reach.

TePe Compact Tuft™ has a firm rounded tuft suitable for **overdentures**.

Interdental brushes

Toothbrushing cleans the buccal and lingual of the implant surfaces. To reach the sides of the implants, an interdental brush is often recommended. TePe's interdental brushes are available in different sizes, texture and design to fit individual needs. The handle gives a comfortable grip and all sizes have a **plastic coated wire for safe use** when cleaning around implants.

TePe Angle™ also available.



TO ORDER FOR YOUR PRACTICE

	Brush	Code	Packaging Box
1	Implant Care	472-368	14 each
2	Implant Brush	442-100	25 each
3	Compact Tuft	432-100	25 each

Welcoming New Professionals

Go to: www.tepecanada.com Select: **Create an Account**

For more information please call **TePe Canada 888-789-9928** or Email: info@tepecanada.com

