



As Featured by Gilberto Debelian, DMD, PhD in his Master Clinical Workshops and Lectures


## ENDODONTIC FILE FORCEPS

- Reduces incidence of joint fatigue when using K-files
- Quick connect/disconnect
- Greater ergonomics
- Two distinct models to choose from:


**PCF-N-90AHF/L**



**PCF-N-75CHF/L**



The sharp prongs may impale the parachute holes facilitating 270°s of positioning in the X, Y and Z axis for easy approach to all canals.



Ideal for use with mirrors or microscopes. For obturation they are ideal for use with Gutta-Percha and resin based carriers.

## STEIGLITZ FORCEPS

North/South angulations address the bucco-lingual chamber orientation of anteriors and bicuspsids.

Available in 45°, 75° and 90°

- With carbide inserts (extra grip)
- With diamond dust (extra thin)



East/West angulations address the mesio-distal chamber orientation of molars.

Available in 45°, 75° and 90°

- With carbide inserts (extra grip)
- With diamond dust (extra thin)

