



As featured by Richard Rubinstein, DDS, MS in the Master Clinician Series

Suture Cutting Needle Holder CE2-731-10RL



Corn Forceps PCF-N-CORN

Corn forceps are used to help guide the precise tissue penetration by the suture needle.



FOR REMOVING EASILY ACCESSIBLE SUTURES

Scissors / Forceps combination N-103F

Scissors/forceps combination for removing suture with one hand/one quick and easy motion...

*Engage the suture on the right side of the knot, close, hold firmly and lift.
The suture is then removed quickly and efficiently.*



MPF-N-4CXF

EH Micropoint scissors are used for removing buried sutures atraumatically.



Gingival Retraction Handle GRH

Specially coated retraction device for atraumatic tissue manipulation.

