



Implant Removal Screws

Universal Implant Removal Screws available for implant removal systems:

- Hubermed Implant Removal Kit Gen II
- Neobiotech Fixture Removal FR Kit
- Zimmer Biomet Implant Removal IR Kit



Screw size	Product Name
M 1.4	IRS 14
M 1.54*	IRS 154
M 1.6	IRS 16
M 1.8	IRS 18
No 1-72 UNF	IRS 172
M 2.0	IRS 20
M 2.5	IRS 25

*M 1.54 exclusively available at Hubermed

Note: Screws in Implant CLOCKWISE

DISPOSABLE: NOT TO BE REUSED TO AVOID FRACTURE WHEN EXPLANTING

Removal Screws are implant specific.
Refer to compatibility list available at www.hubermed.com

Conveniently packaged in sets of two (2) that provide a backup when needed

To Order
or for more info **HUBERMED** 1-888-789-9928
info@hubermed.com