



# IMPLANT REMOVAL KIT GEN II

PREVIOUS KIT FOR EXISTING GEN I USERS

Complete Tool Kit  
Designed To Remove  
Multiple Implant Systems




CONVERSION KIT AVAILABLE  
FOR NEW GEN II



## A new patent-pending WAVE TECHNOLOGY

- The two wave-like shape parts spread out to exert greater downward force due to a new and improved design that reduces rotational slippage to allow for greater control during the implant removal procedure.
- The advantage of 2 wave-like shape parts allows faster constant separation of components after the implant is removed
- Eliminates the use of trephine drills and thereby minimizes bone trauma



### Complete kit includes:

- 12 screws, 2 implant remover bodies,
- 4 implant remover tips, 2 implant remover adapter,
- 3 implant remover drivers, 1 torque wrench

### INCLUDES:

1. Tray and Lid
2. Implant remover body X 2
3. Implant remover tip X 4
4. Implant remover adapter X 2

Promotional offer ends June 30, 2017

**HUBERMED**  
1-888-789-9928  
www.hubermed.com

| Implant Remover Screw | Screw size                   |
|-----------------------|------------------------------|
|                       | M 1.4                        |
|                       | M 1.54*(not included in kit) |
|                       | M 1.6                        |
|                       | M 1.8                        |
|                       | No 1-72 UNF                  |
|                       | M 2.0                        |
|                       | M 2.5                        |

\*M 1.54 exclusively available at Hubermed

