


Compatible with Zest Locator™



For major Implant Systems:

- Astra Tech
- Straumann
- Branemark
- Megagen AnyRidge and more...
- Biomet 3i
- Nobel

Ease of handling

Pre-loaded abutment with disposable friction grip abutment placement carrier included



Dual retention metal housing

119 \$

* Replacement Male Cap Set **Included** with each abutment purchased

Replacement Male Cap Set

Metal Housing White Spacer 2.4 lbs 1.2 lbs Extra Light Angled

included in a capsule with each abutment purchase
* \$35 Also sold separately (Order code: CDPH)

15° angulated abutment

165 \$

* Replacement Male Cap Set **Included** with each abutment purchased

Restorative Components and ALL IN ONE "Magic Tool" one piece construction



Insertion tip used for nylon placement Triangular design used to place or remove abutments Tip used for nylon removal



KERATOR Abutment Analog	DKA3854	Sold in packs of 4	\$60
KERATOR Coping	DKI4845	Sold in packs of 4	\$60
KERATOR Magic Tool	KMT002		\$99
KERATOR Torque Tip	KMD519		\$69
KERATOR Hand Torque Housing	KMH001		\$70

Replacement Nylons - pack of 4 **\$19**

	4 lbs	CPW (white)
	Angled (for lab) 3 lbs	CPG (green)
	2.4 lbs	CPP (pink)
	1.5 lbs Light Angled	CPO (orange)
	1.2 lbs	CPB (blue)
	0.5 lbs Extra Light Angled	CPR (red)
	0 lb NEW ZERO retention	CPY (yellow)
	Spacer NEW pack of 4	DS (white) \$39

KERATOR

Overdenture Attachment Systems

CODE	Manufacturer	Implant Description	Diameter	T.Height
PLEASE CALL FOR "BAR TYPE"				
AS401-406	ASTRA TECH	4.0(ST)	4.0mm	1mm-6mm
AS4AN	ASTRA TECH	4.0(ST)	4.0mm	15degree
AS3508-AS356	ASTRA TECH	3.5	3.5mm	0.8mm-6mm
AS35AN	ASTRA TECH	3.5	3.5mm	15degree
AT401-AT406	ASTRA TECH	1.6 Thread Aqua 3.5, 4.0	4.0mm	1mm-6m
AT4AN	ASTRA TECH	1.6 Thread Aqua 3.5, 4.0	4.0mm	15degree
AS501-AS506	ASTRA TECH	4.5, 5.0(ST) / Lilac 5.0	5.0mm	1mm-6mm
AS5AN	ASTRA TECH	4.5, 5.0(ST) / Lilac 5.0	5.0mm	15degree
BHM351-BHM356	BIO-HORIZON	MAESTRO 3.5 (yellow)	3.5mm	1mm-6mm
BHM35AN	BIO-HORIZON	MAESTRO 3.5 (yellow)	3.5mm	15degree
BHM401-406	BIO-HORIZON	MAESTRO 4.0 (green)	4.0mm	1mm-6mm
BHM4AN	BIO-HORIZON	MAESTRO 4.0 (green)	4.0mm	15degree
BHM501-506	BIO-HORIZON	MAESTRO 5.0 (blue)	5.0mm	1mm-6mm
BHM50AN	BIO-HORIZON	MAESTRO 5.0 (blue)	5.0mm	15degree
BHP571-576	BIO-HORIZON	PRODIGY 5.7 (blue)	5.7mm	1mm-6mm
BHP57AN	BIO-HORIZON	PRODIGY 5.7 (blue)	5.7mm	15degree
SP371-SP376	IMPLANT DIRECT	3.7 SCREWPLANT	3.7mm	1mm-6mm
SP37AN	IMPLANT DIRECT	3.7 SCREWPLANT	3.7mm	15degree
SP471-SP476	IMPLANT DIRECT	4.7 SCREWPLANT	4.7mm	1mm-6mm
SP47AN	IMPLANT DIRECT	4.7 SCREWPLANT	4.7mm	15degree
SP571-SP576	IMPLANT DIRECT	5.7 SCREWPLANT	5.7mm	1mm-6mm
SP57AN	IMPLANT DIRECT	5.7 SCREWPLANT	5.7mm	15degree
B401-406	Biomet (3l)	4.0 (HEX)	4.0mm	1mm-6mm
B4AN	Biomet (3l)	4.0 (HEX)	4.0mm	15degree
IN501-506	Biomet (3l)	4.0 (HEX)	4.0mm	1mm-6mm
IN50AN	Biomet (3l)	5.0 (HEX)	5.0mm	15degree
CE400-406	Biomet (3l)	Certain 4.1	4.0mm	0mm-6mm
CE4AN	Biomet (3l)	Certain 4.1	4.0mm	15degree
T400-406	ITI STRAUMANN	Soft tissue level 4.8	4.0mm	0mm-6mm
T4AN	ITI STRAUMANN	Soft tissue level 4.8	4.0mm	15degree
CF331-336	ITI STRAUMANN	Bone level 3.3	3.3mm	1mm-6mm
CF33AN	ITI STRAUMANN	Bone level 3.3	3.3mm	15degree
CF401-406	ITI STRAUMANN	Bone level 4.0	4.0mm	1mm-6mm
CF4AN	ITI STRAUMANN	Bone level 4.0	4.0mm	15degree
AR401-406	MEGAGEN	AnyRidge	All	1mm-6mm
AR4AN	MEGAGEN	AnyRidge	All	15degree

CODE	Manufacturer	Implant Description	Diameter	T.Height
PLEASE CALL FOR "BAR TYPE"				
B331-336	NOBEL BIOCARE	BRANEMARK NP	3.3 mm	1mm-6mm
B33AN	NOBEL BIOCARE	BRANEMARK NP	3.3 mm	15degree
B401-406	NOBEL BIOCARE	BRANEMARK RP	4.0 mm	1mm-6mm
B4AN	NOBEL BIOCARE	BRANEMARK RP	4.0 mm	15degree
B501-506	NOBEL BIOCARE	BRANEMARK WP	4.0 mm	1mm-6mm
B5AN	NOBEL BIOCARE	BRANEMARK WP	4.0 mm	15degree
SH3251-3256	NOBEL BIOCARE	HEX-LOCK (HL) 3.25 (HEX)	3.25mm	1mm-6mm
SH325AN	NOBEL BIOCARE	HEX-LOCK (HL) 3.25 (HEX)	3.25mm	15degree
SR351-6mm	NOBEL BIOCARE	REPLACE 3.5 (HEX)	3.5mm	1mm-6mm
SR35AN	NOBEL BIOCARE	REPLACE 3.5 (HEX)	3.5mm	15degree
SR501-06	NOBEL BIOCARE	REPLACE 5.0 (HEX)	5.0mm	1mm-6mm
SR50AN	NOBEL BIOCARE	REPLACE 5.0 (HEX)	5.0mm	15degree
CC400-06	NOBEL BIOCARE	STERIOSS 3.8	4.0mm	0mm-6mm
CC40AN	NOBEL BIOCARE	STERIOSS 3.8	4.0mm	15degree
SRS350-356	NOBEL BIOCARE	REPLACE SELECT 3.5	3.5mm	0mm-6mm
SRS35AN	NOBEL BIOCARE	REPLACE SELECT 3.5	3.5mm	15degree
SRS430-436	NOBEL BIOCARE	REPLACE SELECT 4.3	4.3mm	0mm-6mm
SRS43AN	NOBEL BIOCARE	REPLACE SELECT 4.3	4.3mm	15degree
SRS501-506	NOBEL BIOCARE	REPLACE SELECT 5.0	5.0mm	0mm-6mm
SRS50AN	NOBEL BIOCARE	REPLACE SELECT 5.0	5.0mm	15degree
NA350-356	NOBEL BIOCARE	NOBELACTIVE 3.5	3.5mm	0mm-6mm
NA35AN	NOBEL BIOCARE	NOBELACTIVE 3.5	3.5mm	15degree
NA450-456	NOBEL BIOCARE	NOBELACTIVE 4.3, 5.0	4.0mm	0mm-6mm
NA45AN	NOBEL BIOCARE	NOBELACTIVE 4.3, 5.0	4.0mm	15degree
CO401-406	ZIMMER - CALCITEK	OMNILOC 4.0	4.0mm	1mm-6mm
CO40AN	ZIMMER - CALCITEK	OMNILOC 4.0	4.0mm	15degree
CI3250-256	ZIMMER - CALCITEK	OMNILOC 3.25	3.25mm	0mm-6mm
CS325AN	ZIMMER - CALCITEK	OMNILOC 3.25	3.25mm	15degree
PA350-356	ZIMMER - PARAGON	SCREW-VENT 3.3, 3.7, 4.1	3.5mm	0mm-6mm
PA35AN	ZIMMER - PARAGON	SCREW-VENT 3.3, 3.7, 4.1	3.5mm	15degree
PA450-456	ZIMMER - PARAGON	SCREW-VENT 4.5	4.5mm	0mm-6mm
PA45AN	ZIMMER - PARAGON	SCREW-VENT 4.5	4.5mm	15degree
PA550-556	ZIMMER - PARAGON	SCREW-VENT 5.5 6.0	6.0 mm	0mm-6mm
PA55AN	ZIMMER - PARAGON	SCREW-VENT 5.5 6.0	6.0 mm	15degree
AK401-406	Dentsply Friadent	Ankylos 3.5,4,5,5,5	4.0 mm	1mm-6mm
AK4AN	Dentsply Friadent	Ankylos 3.5,4,5,5,6	4.0 mm	15degree
AKCX401-406	Dentsply Friadent	AnkylosCX	4.0 mm	1mm-6mm
FX341-346	Dentsply Friadent	Frialit -2 &XIVE 3.4	3.4mm	1mm-6mm
FX34AN	Dentsply Friadent	Frialit -2 &XIVE 3.4	3.4 mm	15degree
FX381-386	Dentsply Friadent	Frialit -2 &XIVE 3.8	3.8 mm	1mm-6mm
FX38AN	Dentsply Friadent	Frialit -2 &XIVE 3.8	3.8 mm	15degree
FX451-456	Dentsply Friadent	Frialit -2 &XIVE 4.5	4.5 mm	1mm-6mm
FX51-556	Dentsply Friadent	Frialit -2 &XIVE 5.5	5.5 mm	1mm-6mm

* Others available

Replacement Male

KeraTor Metal Housing CDPH (octa design)

- 4 lbs CPW (white)
- Angled (for lab) 3 lbs CPG (green)
- 2.4 lbs CPP (pink)
- 1.5 lbs Light Angled CPO (orange)
- 1.2 lbs CPB (blue)
- .5 lbs Extra Light Angled CPR (red)
- 0 lb NEW ZERO retention CPY (yellow)
- Spacer NEW (pack of 4) DS (white)


Replacement Male Cap Set



* Also sold separately and included with each abutment purchase

Restorative Components and all purpose "Magic Tool"

- KERATOR Abutment Analog DKA3854
- KERATOR Coping DKI4845
- KERATOR Magic Tool KMT002
- KERATOR Torque Tip KMD519
- KERATOR Hand Torque Housing KMH001



* Replacement Male Cap Set **Included** with each abutment purchased

To learn more about this product, please visit us online at www.hubermed.com or to speak to a sales representative, call I (888) 789.9928

Ordering information

Catalog Number Description

SEE ABOVE

* May not be available in all countries

HUBERMED is the exclusive Canadian distributor for KeraTor Overdenture Attachment Systems

* Registered trademarks: Biohorizon, Astra, Ankylos Dentsply, Zimmer® Dental, Nobel Biocare™, Straumann® Dental Implant System, Implant Direct™, Biomet 3®, Megagen AnyRidge®



HUBERMED inc.

Tel: I (888) 789.9928

E-mail: info@hubermed.com

www.hubermed.com