

CUTTING EDGE TECHNOLOGY in **STEIGLITZ STYLE ENDODONTIC FORCEPS**

For retrieving separated files and silver points. Will not splay, slip off or scissor when only a minimal purchase point is available.
All forceps are micro-tipped.



- North/South angulations address the bucco-lingual chamber orientation of anteriors and bicuspsids.
- East/West angulations address the mesio-distal chamber orientation of molars.

With serrated carbide inserts

N/S Forceps with thumb lock

PCF-N-45SL/M

PCF-N-75SL/M

PCF-N-90SL/M

E/W Forceps with thumb lock

PCF-N-45SPL/M

PCF-N-75SPL/M

PCF-N-90SPL/M



OPEN

LOCKED

PROFILE

Extra thin for deeper access

N/S Micro Diamond Dusted Forceps with thumb lock

PCF-D-N-45SL/M

PCF-D-N-75SL/M

PCF-D-N-90SL/M

E/W Micro Diamond Dusted Forceps with thumb lock

PCF-D-N-45SPL/M

PCF-D-N-75SPL/M

PCF-D-N-90SPL/M

