

ENDODONTIC FILE FORCEPS

- Reduces incidence of joint fatigue when using K-files
- Quick connect/disconnect
- Greater ergonomics
- Two distinct models to choose from

PCF-N-90AHF/L



PCF-N-75CHF/L



PCF-N-90AHF/L




PCF-N-75CHF/L



Forceps may circumnavigate the file for easily accessed canals,

or

the sharp prongs may impale the parachute holes facilitating 270°s of positioning in the X, Y and Z axis for easy access to all canals.



As file forceps they are ideal for use with mirrors or microscopes.
For obturation they are ideal for use with Gutta-Percha and resin based carriers.

***Discovering the problems...
creating the solutions.***